



*Go Anywhere  
Do Anything*



1
















KÜSCHALL ULTRA-LIGHT

The küschall Ultra-Light is a light and compact wheelchair designed to suit an active lifestyle, offering both high mobility and independence.

The rigid front makes the küschall Ultra-Light a very robust wheelchair, giving it a sturdiness and driving performance comparable with a rigid frame wheelchair. With its easy folding system, it allows those with reduced strength and/or hand function to fold the wheelchair without having to apply excessive force.

The küschall Ultra-Light foldable active wheelchair offers the ultimate support when needed. The combination of design, driving performance and manoeuvrability results in an ultra-light wheelchair that delivers on performance.

## TECHNICAL DATA

 Seat width: 280 – 500 mm in increments of 20 mm	 Seat depth: 320 – 500 mm in increments of 20 mm	 Seat height: front: 430 – 530 mm rear: 390 – 500 mm in increments of 10 mm
 Knee-to-heel length: 200 – 510 mm in increments of 10 mm	 Backrest height: 300 – 510 mm in increments of 15 mm	 Backrest angle: 82°/86°/90°/ 94°/98°/102°
 Max. safe slope: 7°	 Reversing width: 920 – 1230 mm	 Total length: 70°: approx. 880 mm 80°: approx. 830 mm
 Total width: SW + 180 mm	 Total width folded: approx. 290 mm	 Max. user weight: 130 kg
 Total weight unloaded: approx. 10.1 kg	 Transport weight w/out rear wheels: approx. 7.2 kg	 Crash Test approved Fulfills the requirements of ISO 7176-19

## LEGEND

- CONFIGURATION:** Product shown has an aluminium standard frame in two-tone frame colour 'dynamite orange' (front) and matt black (cross-brace and rear), carbon clothes-guards, light backrest upholstery, two-piece plastic footplates, 4" Starec front wheels and hub brake rear wheels with Schwalbe tyres and aluminium grey handrims.
- FRAME:** Two front frames are available on the küschall Ultra-Light: active 70° or dynamic 80°. Both are also available in a tapered version. The latter option reduces the front width for better manoeuvrability and increased leg support.
- BACKREST:** The küschall Ultra-Light has an optional foldable backrest. Combined with the foldable frame, this makes it easy to handle, move and transport. Measuring just 290 mm wide when folded, the küschall Ultra-Light fits into the smallest of storage spaces.
- HUB BRAKES:** With integrated brakes on the hubs, 24 silver spokes and a black rim, the new hub brake wheels eliminates the need standard brakes and clamps to be fitted. Improved sleek frame design and significant weight saving.
- SEAT BAG:** The küschall Ultra-Light comes with an optional seat bag that can be attached with fixation clips. Placed just in front and below the seat, this bag is situated in an area that is very accessible.
- SWISS-MADE:** All küschall wheelchairs are manually assembled in our factory in Switzerland and have a 'Swiss Made' label attached to the frame.

**United Kingdom:** Invacare Ltd. | Pencoed Technology Park, Pencoed, Bridgend CF35 5AQ | Phone +44 (0) 1656 776222 | Fax +44 (0) 1656 776220 | [www.invacare.co.uk](http://www.invacare.co.uk)  
**Ireland:** Invacare Ltd. | Unit 5, Seatown Business Campus | Seatown Road | Swords, Co. Dublin | Phone +353 (0)1 810 7084 | Fax +353 (0)1 810 7085 | [www.invacare.ie](http://www.invacare.ie)

küschall is a registered trademark. Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Küschall AG reserves the right to modify existing models at any time. Any use of this publication or of parts thereof as well as any reproduction of images must have the written consent of Küschall AG. Patents pending  
 Ultra-Light\_UK\_2014\_1486602 Copyright © 10/2014, Küschall AG, Switzerland - All rights reserved



2



3



4



5



6



KÜSCHALL ULTRA-LIGHT