

rea |  **Azalea Tall**



### **No limits on comfort for taller users**

*Rea **Azalea Tall*** is specially designed to meet the needs of tall users who require a Tilt in Space wheelchair with a longer seat support. Adapted from the *Rea **Azalea***, the *Rea **Azalea Tall*** boasts all the advantages of a reliable, tilting wheelchair and offers a unique weight-shifting mechanism. The *Rea **Azalea Tall***: at the top of its class for comfort and stability.



**Yes, you can.®**



**Features and Options**

**Real compatibility**

The *Rea Azalea Tall* has all the features of the *Rea Azalea* such as Tilt in Space weight-shifting for a really stable tilt. All accessories such as rear wheels for self propelling or smaller wheels for transport and a different seat system also fit the *Rea Azalea Tall*, which is in fact, simply an adaptor kit for the standard *Rea Azalea* model.

**Easy maneuverability**

Attendants will appreciate the adjustable push bar, which offers good ergonomic grip and can be adopted to fit different personal length's of the attendant.



**Longer seat plate**

An extra 5 cm of seat depth is designed to support longer legs.



**Chassis length**

An extra 5 cm in length enables greater forward stability.



**Seat height**

An extra 5 cm of height allows for longer leg-rests.

**Technical data**

<i>Rea Azalea Tall</i>	39 - 59 cm	48 - 55 cm	50 cm	50 - 80 cm (without seat cushion)	24 - 36 cm	33 - 50 cm	-1° - 20° (22/24") -1° - 25° (16")
<i>Rea Azalea Tall</i>	0°-30°	SW + 22 cm (SW + 23 cm)	100 - 130 cm	95 - 107 cm	34 kg (SW 44 cm)		
<i>Rea Azalea Tall</i>	Max 135 kg	Transport weight 18,7 kg	Crash-tested ISO 7176-19	<p>The <i>Rea Azalea Tall</i> complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information regarding vehicular transportation and for more pre-sale information for this wheelchair, please consult the following pages on our website: <a href="http://www.invacare-rea.com">www.invacare-rea.com</a> under Rea wheelchairs/Technical downloads.</p>			

Invacare International Sarl  
Route de Cité-Ouest 2 - 1196 Gland Switzerland  
Tel: +41 (0)22 354 60 10  
Fax: +41 (0)22 354 60 11  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

© 2010 Invacare® International Sarl  
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Rea Azalea Tall - EU - 08/2010

**Frame colour**

**Upholstery**



Pearl Grey



Grey Dartex TR23



Grey plush TR 32

Rea, Rea design and DSS (Dual Stability System) design are registered trademarks of Invacare® International.



**Yes, you can.®**

