



Invacare[®] | ***Propad***[™] ***Low Profile Cushion***

Low-profile castellated cut-foam for comfort and support

Trusted Solutions,
Passionate People



The *Invacare Propad Low Profile cushion* a lower cushion, rather than a standard height of the Propad Original cushion. Some individuals prefer a low profile cushion, finding it more compatible with their daily living activities.

Features and Options



Cover helps to prevent a build up of shear and friction forces

Cover is two-way stretch, water-resistant and vapour permeable. Cover features an anti-slip base.

Available in a range of sizes

The *Invacare Propad Low Profile Cushion* is available in 102 sizes, ranging from 250 x 250 mm – 600 x 550 mm.

Full length flap conceals the zip

The full length flap concealing the zip reduces the potential for fluid ingress into the cushion and protects the foam core. The cover can be removed and laundered.

Castellated Foam Cushion





Each castellation delivers enhanced support and pressure area care to the individual.

Reduces shear and friction forces

The castellations in the *Invacare Propad Low Profile cushion* reduce shear and friction forces on tissue.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see the [Invacare website](#).

Size range	Max. user weight	Max. cleaning temperature	Warranty
			
A range of 102 sizes	76 kg (12 stone)	80°	1 year

Propad Low Profile Cushion

Manufactured to European Standards:

Interior Foam: BS EN 597-1 Cig, BS EN597-2 Match, BS 5852 (Crib 5)

Cover: BS 7175 (Crib 5)

Colours

Blue PU

Black Towelling



Tel: +44 1 656 776 222
 Fax: +44 1656 776 220
 E-mail: uk@invacare.com
 Sales Order E-mail: ordersuk@invacare.com
www.invacare.co.uk

Invacare Limited
 Pencoed Technology Park
 Pencoed - Bridgend CF35
 5AQ - United Kingdom

Tel. +353 1 8107084
 Fax +353 1 8107085
 Email: ireland@invacare.com
www.invacare.ie

Invacare Ireland Ltd
 Unit 5 - Seatown
 Business Campus
 Seatown Road
 SWORDS - County
 Dublin - Ireland



Yes, you can.®

