




### Black Smooth




### Black wipedown



### Black spacer fabric



### Baltic Blue



#### Smooth Cover

*(Also known as Ingress Resistant Cover)*

- A black fabric material
- Ingress-resistant treated fabric
- Multi stretch
- Vapour-permeable
- Reduces potential for perspiration and irritation of the skin
- Facilitates easy transfers

#### Available with:

Flo-tech Contour, Flo-tech Contour Lo-back, Flo-tech Contour Visco, Flo-tech Image, Flo-tech Lite, Flo-tech Lite Visco, Flo-tech Plus, Flo-tech Solution, Flo-tech Solution Xtra, Flo-tech Deep Solution, Essential Basic, Essential Visco, Propad Original, Propad Premier, Propad Profile

#### Wipedown Cover

- A black fabric outer surface
- Multi-stretch
- Water resistant, easy to maintain wipedown surface
- Vapour-permeable
- Facilitates easy transfer
- Suitable for clients with continence issues

#### Available with:

Flo-tech Contour, Flo-tech Contour Lo-back, Flo-tech Contour Visco, Flo-tech Image, Flo-tech Lite, Flo-tech Lite Visco, Flo-tech Plus, Flo-tech Solution, Flo-tech Solution Xtra, Flo-tech Deep Solution, Essential Basic, Essential Visco

#### Spacer Cover

- Black cover
- Allows air to circulate
- Reduces the risk of moisture/heat build up
- Multi-stretch
- May prevent maceration

#### Available with:

Aerotech Uno, Aerotech 4, Flo-tech Contour, Flo-tech Contour Lo-back, Flo-tech Contour Visco, Flo-tech Image, Flo-tech Lite, Flo-tech Lite Visco, Flo-tech Plus, Flo-tech Solution, Flo-tech Solution Xtra, Flo-tech Deep Solution

#### Two way Stretch PU

- Available in Baltic Blue
- Water resistant
- Wipedown surface
- Certain cushion types have limited colour options

#### Available with:

Flo-tech Contour, Flo-tech Contour Visco, Flo-tech Deep Solution, Flo-tech Image, Flo-tech Lite, Flo-tech Lite Visco, Flo-tech Plus, Flo-tech Solution, Flo-tech Solution Xtra



Max. cleaning temperature



80°

Warranty



3 years

Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [ireland@invacare.com](mailto:ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

Invacare Ireland Ltd  
Unit 5 - Seatown  
Business Campus  
Seatown Road  
SWORDS - County  
Dublin - Ireland